QUICK REFERENCE GUIDE


Walgreens
The Brand America Trusts®
BM-725W

1 BEFORE TAKING A MEASUREMENT

1 Sit quietly in a comfortable environment for about 10 minutes. Try to refrain from talking and moving.

2 Wait 30-45 minutes before taking a measurement if you've just consumed any caffeinated beverages or smoked cigarettes /cigars.

3 For consistency, try to take your blood pressure using the same arm, we suggest the left.


4 Sit with your legs uncrossed and your feet flat on the floor. Do not touch the cuff or monitor while measuring.

5 Wait at least 10 minutes in between measurements.

6 Note: If your arm measures more than 13", please use Walgreens Large Cuff (Model # BC-200W)


2 CUFF PLACEMENT

1 Make sure you are not wearing a thick coat or sweater. If you're wearing a long sleeved shirt,


roll up the sleeve so it's not constricting circulation to your arm.

2 Wrap the cuff around your arm making sure that it's not too tight.

For detailed instructions on how to apply the cuff, refer to p. 16 in the manual.

3 Place your arm flat on a table with your palm facing upwards.

QUICK REFERENCE GUIDE

Walgreens
The Brand America Trusts®
BM-725W

3 TAKING A MEASUREMENT

1 Once the cuff is properly applied and you are seated with your feet flat on the floor you are ready to take a measurement.

2 Press the START/STOP button . When the " " depress knob symbol flashes on the LCD screen, depress the knob to release any trapped air in the cuff. Next, an

upwards arrow will flash indicating you to squeeze the bulb repeatedly to inflate the cuff. *Do not inflate higher than 240mmHg.

3 You will see the numbers quickly change as you inflate the cuff. Pump the bulb until you hear a beep. Do not continue to inflate. The monitor will begin the measurement process as it deflates.

4 Once measurement is complete the monitor will automatically deflate and the " " symbol will appear indicating you to depress the knob to release any remaining trapped air.

5 Your blood pressure and pulse measurements will be displayed on the digital panel. The monitor will automatically power off after about 1 minute.

Walgreens c/o SAI
105 CHALLENGER ROAD, 6FL
RIDGEFIELD PARK, NJ 07660

Toll Free Customer Service
1.800.261.0130

Healthy Living® is a registered trademark used under license by Samsung America, Inc.

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Mejor Salud® están bajo la licencia de Samsung America, Inc.